Ultimaster™ Tansei™

Sirolimus eluting coronary stent system

MASTERING COMPLEXITY.
SIMPLE.

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
Ultimaster™ Tansei™

Sirolimus eluting coronary stent system


	Cardiology portfolio	
	Current portfolio	Recent update
ESM Entry Site Management	Diagnostic catheter Guidewire Introducer sheath Radial band VCD	 Glidesheath Slender Challenge for complex PCI through radial
TLA Therapeutic Lesion Access	Semi-compliant PTCA balloon Non-compliant PTCA balloon PTCA guidewire Guiding catheter Coronary microcatheter OFDI Aspiration catheter	 Accuforce Expandability for tough lesion including calcification  Ryurei Upgrade deliverability for complex PCI
Therapeutic Device	DES	 Ultimaster Tansei Improved deliverability for complex PCI

Ultimaster™ Tansei™

Sirolimus eluting coronary stent system

Ultimaster™ Tansei™

Sirolimus eluting coronary stent system

PUSHING BOUNDARIES

Workhorse practice

Complex PCI

High bleeding risk patients

Complication reduction

Mastering complexity. Simple

Value

Increase success rates from simple to complex cases.
Save time and reduce stress; treat better, treat quickly

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Removal Procedure
Step Action

1. Ensure balloon is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position, withdraw Stent Delivery System.

Note: Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system, the entire system should be removed. See stent/System Removal – Precautions section for removal instructions.

4. Tighten rotating hemostatic valve.
5. Repeat angiography to assess stented area. If an adequate expansion has not been obtained, exchange back to a balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall. Do not dilate the stent beyond the limit tabulated below.

The Ultimaster Tansei Sirolimus Eluting Coronary Stent should not be expanded to a diameter beyond 0.5 mm of the reference vessel diameter.

- 6. Final stent diameter should match reference vessel diameter.

removal instructions.

4. Tighten rotating hemostatic valve.
5. Repeat angiography to assess stented area. If an adequate expansion has not been obtained, exchange to a balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall. Do not dilate the stent beyond the limit tabulated below.

Nominal stent diameter (ID)	Post-Dilatation limits (ID)
Φ2.25mm, 2.5mm, 2.75mm, 3.0mm	4.5mm
Φ3.5mm, 4.0mm	5.5mm

6. Final stent diameter should match reference vessel diameter.

33	DE-RQ2233KSM	DE-RQ2533KSM	DE-RQ2733KSM	DE-RQ3033KSM	DE-RQ3533KSM	DE-RQ4033KSM
38	DE-RQ2238KSM	DE-RQ2538KSM	DE-RQ2738KSM	DE-RQ3038KSM	DE-RQ3538KSM	DE-RQ4038KSM

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

NEWS • Conference MASTER Enough Bleedi

The NEW ENGLAND
JOURNAL of MEDICINE

This article is available to subscribers. [Subscribe now](#). Already have an account? [Sign in](#)

SUBSCRIBE
OR RENEW

ORIGINAL ARTICLE [FREE PREVIEW](#)

Dual Antiplatelet Therapy after PCI in Patients at High Bleeding Risk

Marco Valgimigli, M.D., Ph.D., Enrico Frigoli, M.D., Dik Heg, Ph.D., Jan Tijssen, Ph.D., Peter Juni, M.D., Pascal Vranckx, M.D., Ph.D., Yukio Ozaki, M.D., Ph.D., Marie-Claude Morice, M.D., Bernard Chevalier, M.D., Yoshinobu Onuma, M.D., Ph.D., Stephan Windecker, M.D., Pim A.L. Tonino, M.D., et al., for the MASTER DAPT Investigators*

Search

Clinical Topics Latest In Cardiology Education and Meetings Tools and P

Management of High Bleeding Risk Patients Post Bioresorbable Polymer Coated Stent Implantation With an Abbreviated vs. Standard DAPT Regimen - MASTER DAPT

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